CHLSO External Charger

The CHL50 charger is an important part of the LIFT50 family. It functions as a wall charger when mounted on the wall, but also as a more mobile charger simply placed horizontally on its back on any surface.

Usage

Usage temperature:	+5 °C to 40 °C
Storage temperature:	-10 °C to 50 °C
Relative humidity:	20% to 80% - non-condensing
Atmospheric pressure:	700 to 1060 hPa
Height above sea level:	Max. 3000 meters

Mains supply grid should be limited to the highest prospective short circuit current of 100 A.

Demands to mains supply safety for the application in accordance with IEC 60601-1, ed. 3.2, § 8.11.5

Nominal current draw:	Max. 350 mA
Power consumption (standby):	<0.5W
Power consumption (charging):	Max. 25 W
Approvals:	IEC 60601-1 ANSI/AAMI ES60601-1 CAN/CSA-22.2 No. 60601-1 IEC 60601-1-2 PSE


LED1 - LED2 - LED3	LED state (Not listed = off)	States while charging
	LED 1-3 constantly on	With battery or if no battery mounted on charger: 90-100%
	LED 1+2 constantly on LED 3 flashes slowly	65-90%
	LED 1 constantly on LED 2 flashes slowly	40-65%
	LED 1 flashes slowly	0-40%
	LED 1+2+3 flash slowly	Charging stopped due to low battery temperature, high battery temperature or other error conditions
	No light in LEDs	Charging stopped due to lost communication to battery

Charger mounting

When mounting CHLS0 on wall, please use the two dedicated mounting holes placed in the charger body.

CHLS0 must be mounted by means of 2 screws (not supplied by LINAK):

Screw type:


Diameter 4 mm/style: pan head or truss head.
Length and thread type depend on wall material.

Washer type:

Optional 4 mm washer with max. 10 mm
outer diameter.

The mounting screws for the control box and the charger must be tightened with a maximum torque of 1 Nm.

When placing CHLS0 on table or shelf, 4 rubber studs on the back add stability and friction to the charger/battery combination.

Recommendations - mounting

As long as the charger is mounted correctly, then the CHL50 complies to IPX4.

If the CHL50 is mounted incorrectly, then water will gather in the plug holes resulting in possible non-compliance with IPX4.

CHL50 can be mounted as shown on the pictures:

- 1) Batteries up, cable outlets downwards
- 2) CHL50 lying on the back, no difference if mains cable is up or down
- 3) IPX4 only if hanging upright

Hanging on wall - Vertical

Hanging on wall - Horizontal

Laying flat on any surface

Recommendations

- The charger CHL50 is specifically designed to charge the BAL50 type battery.
- Special care should be taken when mounting the CHL50.
- If the CHL50 is mounted correctly, the CHL50 complies to IPX4.
- If the CHL50 is mounted incorrectly, then water may enter the screw holes resulting in IPX4 non-compliance and cause malfunction and hazardous situations.
- In cleaning situations, the charger must be disconnected from mains.

Warnings

- Make a proper cable installation to avoid short-circuit of cables. Regular inspection must be made for wear and damage. Defective parts must be replaced.
- Ensure that mains cable plug is fully inserted into mains socket.

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