

 TOPRO

Premium rollators



UK

made in Norway



TOPRO Pegasus

Your lightweight premium companion for effortless mobility.



Modern, slender and lightweight carbon rollator with functionality and elegant design. This design combines to guarantee the highest level of quality, safety and stability at a weight of just 6.1 kg*.

Tried and tested features have been further developed:

- User-friendly method for setting the height of the handles using a **numeric scale**.
- **Partly integrated brake cables** and modern reflective panels on the frame.
- Stylish detachable shopping bag, which, with its retracted position provides **increased forward view to improve user safety**.
- Easy handling and manoeuvrability.
- **Reflectors** on the frame, lock wheels and shopping bag provide **360 degree visibility for safety in the dark**.
- **Shock-absorbing tyres and ball bearings** in wheels and forks for a comfortable ride.

Available colours:



Measurements, colour options and other product details are available on pagina 18.



* Weight including shopping bag 6.6 kg.



“TOPRO Pegasus represents a new trend in the development of mobility aids where design, functionality and safety are combined to meet users' expectations.”



made in Norway



Troja 5G

Your stylish companion for more independence in everyday life.



Seamless combination of elegance and function – true to the TOPRO tradition. The TOPRO Troja 5G Black Series impresses with modern, lightweight design and new colour variants. It offers high stability and a safe driving experience with maximum comfort.

BLACK SERIES - *Move in style*

Details make the difference – for greater user-friendliness and comfort:

- Intuitive handle height adjustment with a **numerical scale** and **memory function**
- **Shock-absorbing wheels** provide maximum comfort on various surfaces
- **Modern and innovative shopping bag** with up to 10 kg capacity*
- **Ergonomically shaped handles** (Ergo-Grip)
- Easy handling and manoeuvrability thanks to **ball bearings** in forks and wheels
- Quick-release axle for **easy wheel change**

Dimensions, colours, and further product details can be found on page 14.

Model M & S available in new matt-coated colours:



* The design of the bag differs from the depiction.



Learn more
Troja 5G



*For the precious
moments together.*

TOPRO TROJA 5G -
BLACK SERIES

made in Norway



Olympos ATR (Slim)



Your trusted outdoor companion – for maximum independence.

With extra-large front wheels and off-road tyres, the Olympos ATR Black Series in its new design offers reliable and comfortable support on uneven roads, cobblestones, and unpaved paths. The perfect choice for outdoor activities. Also available in a slimmer version (ATR Slim).

BLACK SERIES -
Move in style

Innovative features for maximum support:

- **Wider frame** for greater stability when used outdoors.
- **Large shock-absorbing tyres** ensure a comfortable ride and improved manoeuvrability on rough terrain.
- **Modern and innovative shopping bag** with up to 15 kg capacity*
- Intuitive handle height adjustment with a **numerical scale** and **memory function**
- Easy driving and manoeuvring thanks to **ball bearings** in forks and wheels

Model M, M Slim & S available in new matt-coated colours:



Dimensions, colours, and further product details can be found on page 15.



The ergonomic handle (**TOPRO ErgoGrip**) makes it easier and safer to stand up from a sitting position.

* The design of the bag differs from the depiction.

Learn more



Olympos ATR



*For unforgettable
moments with
family and friends.*

TOPRO OLYMPOS ATR –
BLACK SERIES

made in Norway



Troja Original



Your reliable companion for an independent life.

The Troja Original Black Series is our quality rollator for beginners, combining modern design, high quality, and pioneering features that have influenced the entire rollator industry. It offers all the functions you need for a safe and active everyday life. Now in a new design.

BLACK SERIES -
Move in style

The functionality is based on TOPRO's continued high quality over the years:

- **Lightweight, only 6.8 kg**, with a high user weight capacity of 165 kg.
- The trend-setting **longitudinal folding** is easy to operate with one hand by pulling the seat strap.
- **Ball bearings** in forks and wheels enable easy driving and manoeuvring.
- **Practical shopping net basket** for direct access to contents.
- Available in three different sizes.
- Easy to handle and operate, with great manoeuvrability, making it a reliable companion in everyday situations.
- **Scale for easier adjustment** of handle height. With improved locking bracket.

Model M, S & X available in new matt-coated colours:



Dimensions, colours, and further product details can be found on page 16.



Learn more



Troja Original



For the joy of movement.

TOPRO TROJA ORIGINAL -
BLACK SERIES

made in Norway



Troja Walker²



Your secure forearm-height companion – for stable walking support.

The lightweight and practical TOPRO Troja Walker² is a stable walking aid with the comfort of a rollator. Its forearm supports assist with upright walking and ensure correct posture.

- Easy-to-use **brake and parking brake** – operated with one brake lever.
- The trend-setting **longitudinal folding** can be operated with one hand by pulling the seat strap. **Easy to transport**, stable standing on the wheels, and space-saving storage.
- Handles can be pushed down into the frame for storage.
- The unique **memory function** makes it easy to find the pre-set handle height.
- **Ball bearings** in forks and wheels enable easy driving and manoeuvring.
- **Edge deflectors** and **tip assist** for safe and easy manoeuvring.
- **Practical shopping net basket** with direct access to items. Load capacity up to 5 kg.
- **Scale for easier handle height adjustment**, with improved safety bracket.
- Available in three different sizes.

Dimensions, colours, and further product details can be found on page 17.

Model M, S & XS available in:



Learn more




Troja Walker²




The TOPRO rollator family for everyday use

Whether it is product features, comfort or design – we have continually improved our rollators and expanded our line-up for over 30 years.

We are convinced our products provide you with an added value – by combining smart functionality and modern design with high quality and safety.



Accessories for added comfort



Crutch holder

Can be fitted on either side.


Troja Original / 5G,
Olympos ATR, grey: Art. No. 815358
 Black Series: *Troja Original / Olympos ATR /*
5G, black: Art. No. 814040
TOPRO Pegasus: Art. No. 815236



Back support with padding

Ergonomically height adjustable, available in two lengths.

TOPRO Pegasus, Troja 5G,
Olympos ATR (Slim):
 Short 67 cm: Art. No. 815270
 Long 75 cm: Art. No. 815272
Troja Original:
 Short 68 cm: Art. No. 815271
 Long 76 cm: Art. No. 815273



Tray

Snaps in securely, easy to remove and clean, carries up to 5 kg.

TOPRO Pegasus,
Troja Original / 5G: Art. No. 814728
Olympos ATR M: Art. No. 815212
Olympos ATR Slim: Art. No. 815363



Anti-slip mat for tray

For safe transport of items.

TOPRO Olympos ATR / Pegasus /
Troja Original /
Walker / 5G: Art.Nr. 814059
Olympos ATR Slim: Art.Nr. 815832



Umbrella with attachment

Keep dry under the umbrella in rainy weather or use it as a parasol.

Black: Art. No. 814068
 Red: Art. No. 814099

TOPRO Pegasus, Troja 5G / Original,
Olympos ATR (Slim)




One-handed brake

Can be fitted on either side for one-handed operation of simultaneous braking.

Art. No. 814026

TOPRO Pegasus, Troja 5G / Original,
Olympos ATR (Slim)



LED lamp with clip on and USB plug for safety in the dark.
 Chargeable.

Art. No. 815260

TOPRO Pegasus, Troja 5G / Original,
Olympos ATR (Slim)



Rear bag with zipper

Can be opened with one hand, close to the body, ideal for your valuables. Carries up to 3 kg.

Art. No.

814046



TOPRO Pegasus, Troja 5G / Original,
Olympos ATR (Slim)

You can find further accessories on our website.



Technical details



Model	Troja 5G	Troja 5G
Size	M (medium)	S (small)
Recommended user height	150–200 cm	135–170 cm
Max. user weight	150 kg	125 kg
Max. load shopping bag / mesh basket	10 kg / 15 l	10 kg / 15 l
Handle height	78–100 cm	69–86 cm
Handles	Ergonomic, multiple grip positions (TOPRO ErgoGrip)	
Brake lever reach	7 cm	7 cm
Distance between handle centres	47 cm	47 cm
Seat width	45,5 cm	45,5 cm
Seat height	63 cm	54 cm
Turning circle	84 cm	84 cm
Max. height	103,5 cm	89,5 cm
Max. width	60 cm	60 cm
Max. length	65 cm	65 cm
Wheel diameter front/rear	20 cm	20 cm
Wheel width front/rear	3,6 cm	3,6 cm
Folded length	65 cm	65 cm
Folded width	24 cm	24 cm
Folded height	79,5 cm	73 cm
Transport dimensions folded (L × W × H)	65 × 24 × 73 cm*	65 × 24 × 66 cm*
Product weight (without wheels ¹ and shopping bag ²)	5,6 kg ¹⁺²	5,4 kg ¹⁺²
Product weight (with wheels and shopping bag)	7,8 kg	7,6 kg
Area of use	Indoor & outdoor use	
Art. No.	814600 + Var. No.	814601 + Var. No.
Scan the QR code for technical illustrations on our website		

Colour options



*(excluding wheels and shopping bag)

Technical details



Model	Olymos ATR Slim	Olymos ATR	Olymos ATR
Size	Slim M (medium)	M (medium)	S (small)
Recommended user height	150–195 cm	150–195 cm	140–170 cm
Max. user weight	150 kg	150 kg	150 kg
Max. load shopping bag	15 kg / 10 l	15 kg / 12 l	15 kg / 10 l
Handle height	76–95,5 cm	76–95,5 cm	69–86,5 cm
Handles	Ergonomic, multiple grip positions (TOPRO ErgoGrip)		
Brake lever reach	7 cm	7 cm	7 cm
Distance between handle centres	44 cm	49 cm	50 cm
Seat width	40 cm	47 cm	47 cm
Seat height	60 cm	60 cm	54 cm
Turning circle	97 cm	97 cm	97 cm
Max. height	99 cm	99 cm	90 cm
Max. width	62 cm	68 cm	68 cm
Max. length	77 cm	77 cm	77 cm
Wheel diameter front/rear	28 cm / 23 cm	28 cm / 23 cm	28 cm / 23 cm
Wheel width front/rear	4,2 cm	4,2 cm	4,2 cm
Folded length	77 cm	77 cm	77 cm
Folded width	27 cm	27 cm	27 cm
Folded height	82 cm	87,5 cm	84 cm
Transport dimensions folded (L × W × H)	77 × 27 × 82 cm	77 × 27 × 87,5 cm	77 × 27 × 84 cm
Product weight (without wheels ¹ and shopping bag ²)	–	–	–
Product weight (with wheels and shopping bag)	8,3 kg	8,9 kg	8,7 kg
Area of use	Outdoor use		
Art. No.	814320 + Var. No.	814300 + Var. No.	814305 + Var. No.
Scan the QR code for technical illustrations on our website			

Colour options



The differences in dimensions between the Standard Medium and Slim Medium sizes.

Technical details



Model	Troja Original	Troja Original	Troja Original
Size	M (medium)	S (small)	X (extra small)
Recommended user height	150–200 cm	130–170 cm	120–160 cm
Max. user weight	165 kg	165 kg	130 kg
Max. load shopping net	10 kg / 20 l	10 kg / 20 l	5 kg / 11 l
Handle height	78–100 cm	69–86 cm	64–79 cm
Handles	Anatomical		
Brake lever reach	6,5 cm	6,5 cm	6,5 cm
Distance between handle centres	47 cm	47 cm	37 cm
Seat width	45,5 cm	45,5 cm	36 cm
Seat height	63 cm	54 cm	51 cm
Turning circle	84 cm	84 cm	79 cm
Max. height	101 cm	87 cm	80 cm
Max. width	60 cm	60 cm	51 cm
Max. length	65 cm	65 cm	65 cm
Wheel diameter front/rear	20 cm	20 cm	20 cm
Wheel width front/rear	3,6 cm	3,6 cm	3,6 cm
Folded length	65 cm	65 cm	65 cm
Folded width	24 cm	24 cm	24 cm
Folded height	80 cm	74 cm	66 cm
Transport dimensions folded (L × W × H)	65 × 24 × 80 cm	65 × 24 × 74 cm	65 × 24 × 65 cm
Product weight (without wheels ¹ and shopping bag ²)	–	–	–
Product weight (with wheels and shopping bag)	6,8 kg	6,6 kg	6,1 kg
Area of use	Indoor & outdoor use		
Art. No.	815400 + Var. No.	815410 + Var. No.	815420 + Var. No.
Scan the QR code for technical illustrations on our website			

Colour options



Technical details



Model	Troja Walker ²	Troja Walker ²	Troja Walker ²
Size	M (medium)	S (small)	X (extra small)
Recommended user height	150–185 cm	130–170 cm	125–160 cm
Max. user weight	135 kg	125 kg	80 kg
Max. load shopping net	5 kg / 20 l	5 kg / 20 l	5 kg / 11 l
Height of forearm support points	99–114 cm	85–100 cm	81–96 cm
Handles	Round, Ø 3,5 cm		
Brake lever reach	7,4 cm	7,4 cm	7,4 cm
Distance between handle centres	45 cm	45 cm	33 cm
Seat width	45,5 cm	45,5 cm	36 cm
Seat height	63 cm	54 cm	51 cm
Turning circle	84 cm	84 cm	79 cm
Max. height	124 cm	110 cm	106 cm
Max. width	61 cm	61 cm	51 cm
Max. length	65 cm	65 cm	65 cm
Wheel diameter front/rear	20 cm	20 cm	20 cm
Wheel width front/rear	3,6 cm	3,6 cm	3,6 cm
Folded length	65 cm	65 cm	65 cm
Folded width	33 cm	33 cm	33 cm
Folded height	102 cm	88 cm	85 cm
Transport dimensions folded (L × W × H)	65 × 33 × 102 cm	65 × 33 × 88 cm	65 × 33 × 85 cm
Product weight (without wheels¹ and shopping bag²)	–	–	–
Product weight (with wheels and shopping bag)	8,9 kg	8,6 kg	8,3 kg
Area of use	Indoor use		
Art. No.	815430	815440	815450
Scan the QR code for technical illustrations on our website			

Colour options



Technical details



Model	Pegasus
Size	M (medium)
Recommended user height	150–200 cm
Max. user weight	130 kg
Max. load shopping bag	10 kg / 15 l
Handle height	78–97,5 cm
Handles	Ergonomic, multiple grip positions (TOPRO ErgoGrip)
Brake lever reach	7 cm
Distance between handle centres	47 cm
Seat width	43,5 cm
Seat height	63 cm
Turning circle	84 cm
Max. height	101 cm
Max. width	64 cm
Max. length	65 cm
Wheel diameter front/rear	20 cm
Wheel width front/rear	3,6 cm
Folded length	65 cm
Folded width	27 cm
Folded height	80,5 cm
Transport dimensions (L × W × H)	65 × 27 × 80,5 cm
Product weight (excl. shopping bag)	6,1 kg
Product weight (incl. shopping bag)	6,6 kg
Area of use	Indoor & outdoor use
Art. No.	815020 + Var. No.

Scan the QR code for technical illustrations on our website



Colour options



Keep on moving – with confidence, style and independence.

Mobility that inspires confidence.

Movement means freedom – and how we feel whilst moving makes all the difference.

A rollator from TOPRO is more than just a mobility aid: it's your trusted daily companion for greater comfort, safety, and quality of life.

Thoughtfully designed down to the last detail, a premium rollator from TOPRO offers not only exceptional stability but also intuitive handling and smooth rolling comfort – for an independent life.

Quality you can feel. Design that speaks for itself.

As a Norwegian manufacturer with decades of experience, we stand for the highest quality, Scandinavian design, and one clear mission: to make mobility easier, safer, and more stylish.

Our rollators are as individual as the people who use them. Whether compact for everyday life, comfortable off-road, or with a special focus on design – every TOPRO rollator combines thoughtful functionality with long-lasting quality and modern aesthetics.

Thanks to smart features like ergonomic handles, shock-absorbing wheels, user-friendly height adjustment, and an intelligent folding mechanism, our models adapt to your needs with ease – day by day, step by step.


More than function – a feeling of security.

We believe that those who stay mobile stay active. And those who feel supported along the way gain joy and confidence. That's why we place great importance on ensuring you always feel safe and comfortable with a TOPRO rollator.

Our models are low-maintenance, durable, and can be individually customized with a wide range of accessories. Whether in the city, out in nature, or on the go with family and friends – a TOPRO is your reliable companion, wherever life takes you.

– keep on moving

Best regards,
Your TOPRO Mobility Team



Premium products from Norway

Manufacturer:



TOPRO Industri AS
Rambekkvegen 1
2816 Gjøvik
Norway

Phone: (+47) 61 13 46 00

E-mail: info@topromobility.com

Web: www.topromobility.com



TOPRO products are covered by statutory warranty rights in accordance with EU and national law.

In addition, TOPRO offers a manufacturer's warranty of up to 7 years, depending on the product model.

For the TOPRO Taurus forearm walker, the manufacturer's warranty period is 2 years.

Detailed information about the warranty terms can be found in the warranty terms.

In case of warranty claims, please contact your dealer who provided you with the product.

TOPRO walking aids are tested in cooperation with TÜV according to the required standards.

Quality assurance standards

DIN-EN ISO 9001 – Quality

DIN-EN ISO 14001 – Environment

DIN-EN ISO 11199-2 – Rollator

DIN-EN ISO 10993-5 – Biocompatibility

DIN-EN ISO 13485 – Medical devices:

Quality management systems –

Requirements for regulatory purposes

made in Norway



Dealer