


814009

TOPRO

Oxygen basket // Korb für Sauerstoffflasche // Oksygen kurv // Syrgashållare - Assembly instructions - Montageanleitung -

Monteringsanvisning - Model Troja Classic/Original/2G/5G/Neuro M/S, Walker/Walker2 M/S, Olympos/ATR - Art. No. 814009M1-K, 2025-08-14

Content - Inhalt - Innhold



1x 103431


SWL: 5 kg





2x 717029





1x 814009M1



1 Models wit locking wheel/ mit Festellrad/ med Låserad




2





Models with snaparm/ Hebelverstellung/ hurtigjustering:

Olympos ATR



Troja Original



NB!

Make sure that no parts are loose before using the product!

Achtung!

Stellen Sie sicher, dass keine Teile lose sind, bevor Sie das Produkt in Gebrauch nehmen.

NB!

Pass på at ingen deler er løse før du tar produktet i bruk!

OBS!

Kontrollera att inga delar är lösa innan du använder produkten!